

***De novo* sequencing of chimeric spectra using Casanovo**

Yuhui Hong¹, Isha Gokhale¹, Justin Sanders², Melih Yilmaz¹, Gwenneth Straub¹, Wout Bittremieux³, and William Stafford Noble¹²


The authors declare no competing financial interest.

¹*Department of Genome Sciences, University of Washington*


²*Paul G. Allen School of Computer Science and Engineering, University of Washington*

³*Department of Computer Science, University of Antwerp*


DDA Selects and Isolates the Top- k Precursors for Fragmentation



Co-isolated Precursors in DDA Give Rise to Chimeric Spectra




Chimeric Spectra Contain Fragment Ions from Multiple Co-isolated Precursors




Co-isolated precursors

Select top- k
precursors for
fragmentation




Fragment ion peaks from
co-isolated precursors

Every Spectrum Is Potentially Chimeric



Every Spectrum Is Potentially Chimeric — Sensitive Methods Can Reveal Hidden Co-isolated Peptides



Existing *de novo* Methods on DDA Are Not Designed for Chimeric Spectra


Method	Year	Model architecture			Chimeric spectra	Citations
		CNN	Transformer	Others		
DeepNovo (Tran et al. 2017)	2017	✓				564
SMSNet (Karunratanakul et al. 2019)	2019	✓				97
PepNet (Liu et al. 2023)	2023	✓				86
Casanovo (Yilmaz et al. 2022)	2022		✓			117
GraphNovo (Mao et al. 2023)	2023		✓			59
π -HelixNovo (Yang et al. 2024)	2024		✓			46
PointNovo (Qiao et al. 2021)	2021			✓		128
pNovo 3 (Yang et al. 2019)	2019			✓		112

.....

Original table is from [1], where the citations are updated by counting from Google Scholar on May 26, 2026.

[1] Bittremieux, Wout, et al. "Deep learning methods for de novo peptide sequencing." Mass Spectrometry Reviews 45.3 (2026): 507-526.

Casanovo Is a State-of-the-Art Transformer-Based *de novo* Sequencer




[2] Yilmaz, Melih, et al. "Sequence-to-sequence translation from mass spectra to peptides with a transformer model." Nature communications 15.1 (2024): 6427.

Three Key Modifications Enable Casanovo to Directly Predict Chimera

Casanovo-chim


1. Special separator token ":" for dual-peptide output



Three Key Modifications Enable Casanovo to Directly Predict Chimera

Casanovo-chim


1. Special separator token ":" for dual-peptide output
2. Post-processing to correct chimeric charge states




Three Key Modifications Enable Casanovo to Directly Predict Chimera

Casanovo-chim

1. Special separator token ":" for dual-peptide output
2. Post-processing to correct chimeric charge states
3. Permutation-invariant loss function for chimeric output




Permutation-Invariant Loss Function for Chimeric Output



$$\mathcal{L} = \min(\mathcal{H}(p_1:p_2, \hat{y}), \mathcal{H}(p_2:p_1, \hat{y}))$$

mean per-token
cross-entropy loss


We Constructed a Ten-Species Chimeric and Single-Peptide Spectra Dataset Using MSFragger DDA+



Experimental spectra



MSFragger DDA+



Spectra with peptide annotations across ten species


Preprocessing Removes Unreliable Chimeras and Balances Peptide Representation

1. Spectra with 3+ co-isolated precursors excluded.
2. Top 2 spectra retained per unique peptide or chimeric pair (by q-value).
3. Chimeras downsampled to ≤ 10 per unique peptide for balanced representation.

~5M spectra from ~ 1M peptides — 53.3% chimeric




Peptide-Level Partitioning Prevents Data Leakage between Sets





Modified and unmodified peptides are treated as identical.

Precision-Coverage Curves Evaluate Both Accuracy and Completeness of Predictions




Precision-Coverage Curve Can Also be Calculated at Amino Acid Level




Chimeric Predictions Can Exceed Full Coverage — Average precision (AP) Is Evaluated at Coverage ≤ 1



Area under the curve
= average precision



Casanovo-chim Outperforms non-chimeric model at Both Peptide and Amino Acid Levels



Top-Ranked Predictions Match Casanovo's Retention Time Reliability




Bo Wen
Noble Lab


[4] Wen, Bo, et al. "Carafe enables high quality in silico spectral library generation for data-independent acquisition proteomics." Nature Communications 16.1 (2025): 9815.

[5] Wen, Bo, et al. "Carafe2 enables high quality in silico spectral library generation for timsTOF data-independent acquisition proteomics." bioRxiv (2026): 2026-03.


An Iterative Baseline: Casanovo Performs First-Pass Prediction on Observed Spectra



Second-Pass Prediction Is Performed on Residual Peaks After Subtracting Annotated Peaks



Casanovo-Chim Outperforms Two-Pass Casanovo



Takeaway

- We extend Casanovo to directly predict peptides from chimeric spectra.
- At precision 0.8, Casanovo-Chim identifies 17,281 more unique peptides than the non-chimeric baseline — a ~37% lift.
- Casanovo-Chim predictions are validated by consistent retention time.
- Casanovo-Chim outperforms the two-pass Casanovo baseline.

Acknowledgement

University of Washington Noble Lab

William Stafford Noble
Ghulam Murtaza
Ananthan Nambiar
Doudou Yu
Shengqi Hang
Lincoln Harris
Carlo Melendez
Justin Sanders
Bo Wen

Shuhua Xu
Patrick Yu
Xumeng Zhang
Kylie Trousil
Charles Grant
Gwenneth Straub
Priank Dasgupta
Devanshi Desai
Bhavya Krishnan

University of Antwerp Wout Bittremieux

